

Allergen levels in bedroom dust samples found to be un-correlated with baseline symptoms or treatment efficacy of a novel high efficiency particulate air filtration allergen avoidance pillow (PUREZONE™)

Stillerman A MD^{*,†}; Albrecht M[#]; Waldman J[#];

*Allergy & Asthma Specialists, †Clinical Research Institute Inc. of Minneapolis, #Purezone Technologies

ABSTRACT

Introduction: Bedroom dust samples were assessed for allergen levels in a randomized clinical trial involving PUREZONE with the assumption that allergen levels would be representative of overnight subject exposures and, thus, correlate with symptom severity and treatment efficacy.

Methods: In a crossover randomized clinical trial, symptomatic adults (N=35) having perennial allergic rhinoconjunctivitis and a positive skin prick test (dust mite, dog, or cat) received PUREZONE or placebo during 2-week intervals separated by a 1-week washout. Dust samples were collected from bedding and carpeting adjacent to the bed. Allergen levels were quantified via MARIA 5-plex analysis (Indoor Biotechnologies, VA). Correlations in total symptom scores (TSS) [15-point summation of congestion, itch, sneeze, secretion, and eyes] and treatment efficacy were assessed. Approval was obtained from Western IRB and written informed consent obtained from all research subjects.

Results: Significant levels of one or more sensitizing allergens were detected in 43% of subjects' bedrooms, with significant sensitizing exposures totaling 11% for mite group 1 (>2.0 ug/g), 29% for Can f 1 (>2.0 ug/g), and 11% for Fel d 1 (>1.0 ug/g). ANCOVA models found no correlations between baseline TSS vs. sensitizing allergen exposures overnight, upon-waking, during the day, or before bed. Similar models assessing TSS reductions vs. placebo (treatment efficacy) found no correlations between treatment efficacy vs. sensitizing allergen exposures overnight or upon-waking; treatment efficacy during the day or before bed was not significant and, therefore, not assessed.

Conclusions: The results suggest that assessments of subject exposures may not be accurately represented by allergen levels found in dust samples, particularly since significant allergen levels were measured in only a low proportion of bedrooms, yet all subjects were experiencing symptoms and had no other relevant exposures than dog, cat, or dust mite. Furthermore, neither baseline symptoms nor treatment efficacy were dependent on allergen levels, yet PUREZONE was efficacious for reducing overnight symptoms (p<0.001). This lack of treatment efficacy correlation was unexpected given that the benefits of avoidance are normally greater for those subjects with the highest exposures.

INTRODUCTION

Allergen concentrations found in dust samples are often used to classify allergen exposures in a particular environment. In a randomized clinical trial of a personal air filtration allergen avoidance measure (PUREZONE), we sought to determine if baseline symptoms or treatment efficacy were correlated with dust sample allergen concentrations. Dust samples were collected from carpeting or bedding in accordance with NSLAH¹ methods and quantified using MARIA 5-plex² analysis. The complete efficacy data for the randomized clinical trial is published elsewhere³.

PUREZONE ALLERGEN AVOIDANCE THERAPY



- 1 With the touch of a button, bedroom air enters the fan unit.
- 2 Over 99.97% of allergens and pollutants are removed through the charged-fiber HEPA filter.
- 3 Purified air passes through the pillowcase creating a protective, allergen-free breathing zone.

Figure 1. Active and Placebo Treatments (Active Shown Above): The active and placebo treatments were identical devices, with the exception that the placebo device delivered un-filtered air to the breathing zone.

COMPARISON OF AVOIDANCE STUDIES

Device Category	Study	Efficacy vs. Placebo
Personal Air Filtration Systems	Stillerman et al. ³	Allergic Rhinitis: Symptoms & QOL
	Pedroletti et al. ⁶	Asthma: QOL & Exhaled Nitric Oxide
	Zwemer et al. ⁷	Asthma: Symptoms
	Villaveces et al. ⁸	Asthma: Symptoms
	Verrall et al. ⁹	Asthma: Medication Reduction
Multi-Faceted	Moon et al. ¹⁰	Allergic Rhinitis: Symptoms
	Bjornsdottir et al. ¹¹	Allergic Rhinitis: Symptoms & PNIF
	Morgan et al. ¹²	Asthma: Symptoms
Encasement (Isolated)	Terreherst et al. ¹³	(Allergic Rhinitis)
	Frederick et al. ¹⁴	Asthma: Eosinophil
	Woodcock et al. ¹⁵	(Asthma)
	Cloosterman et al. ¹⁶	(Asthma)
	Dharmage et al. ¹⁷	(Asthma)
	Halken et al. ¹⁸	Asthma: Medication Reduction
	Luczynska et al. ¹⁹	(Asthma)
	Reisman et al. ²⁰	Allergic Rhinitis: Symptoms
Room HEPA (Isolated)	Wood et al. ²¹	(Asthma)
	Antoniceilli et al. ²²	(Asthma)

▲ No Benefit ▲ Marginal ▲ Significant
Reported Clinical Benefit

Figure 5. Review of Allergen Avoidance Measure Efficacy: Reported clinical benefits vs. placebo for a select group of studies.

REFERENCES

1. Vejtla PJ, Friedman W, Marker DA, Clickner R, Rogers JW, Viet SM, et al. First National Survey of Lead and Allergens in Housing: survey design and methods for the allergen and endotoxin components. *Environ Health Perspect.* 2002 May; 110(5):527-32.
2. Earle CD, King EM, Tsay A, Pittman K, Saric B, Vailes L, et al. High-throughput fluorescent multiplex array for indoor allergen exposure assessment. *J Allergy Clin Immunol.* 2007 Feb; 119(2):428-433.
3. Stillerman A, Treney N. Efficacy of PUREZONE(TM), a Novel Local High Efficiency Particulate Air (HEPA) Filtration Pillow Device, in Patients With Perennial Allergic Rhinoconjunctivitis (PARC): A Double-Blind, Placebo-Controlled Randomized Clinical Trial (RCT). *Journal of Allergy and Clinical Immunology.* 2009 Feb; 123(2, Supplement 1):S271.
4. Salo PM, Arbes SJ, Crockett PW, Thorne PS, Cohn RD, Zeldin DC. Exposure to multiple indoor allergens in US homes and its relationship to asthma. *J Allergy Clin Immunol.* 2008 Mar; 121(3):678-684.e2.
5. Juniper EF, Rohrbough T, Melzer EO. A questionnaire to measure quality of life in adults with nocturnal allergic rhinoconjunctivitis. *J Allergy Clin Immunol.* 2003 Mar; 111(3):484-90.
6. Pedroletti C, Millinger E, Dahlén B, Söderman P, Zetterström O. Clinical effects of purified air administered to the breathing zone in allergic asthma: A double-blind randomized cross-over trial [Internet]. *Respir Med.* 2009 May 12; [cited 2009 Jul 23] Available from: <http://www.ncbi.nlm.nih.gov/pubmed/19443189>
7. Zwemer RJ, Karibo J. Use of laminar control device as adjunct to standard environmental control measures in symptomatic asthmatic children. *Ann Allergy.* 1973 Jun; 31(6):284-290.
8. Villaveces JW, Rosengen H, Evans J. Use of laminar air flow portable filter in asthmatic children. *Ann Allergy.* 1977 Jun; 38(6):400-404.
9. Verrall B, Mur DC, Wilson WM, Milner R, Johnston M, Dolovich J. Laminar flow air cleaner bed attachment: a controlled trial. *Ann Allergy.* 1988 Aug; 61(2):117-122.
10. Moon JS, Choi SO. Environmental controls in reducing house dust mites and nasal symptoms in patients with allergic rhinitis. *Yonsei Med J.* 1999 Jun; 40(3):238-243.
11. Bjornsdottir US, Jakubindottir S, Runarisdottir V, Juliusson S. The effect of reducing levels of cat allergen (Fel d 1) on clinical symptoms in patients with cat allergy. *Ann Allergy Asthma Immunol.* 2003 Aug; 91(2):189-94.
12. Morgan WJ, Crain EF, Gruchalla RS, O'Connor GT, Kattan M, Evans R, et al. Results of a home-based environmental intervention among urban children with asthma. *N Engl J Med.* 2004 Sep 9; 351(11):1068-80.
13. Terreherst I, Hak E, Oosting AJ, Tempels-Pavlica Z, de Monchy JGR, Brujnzeel-Koomen CAFM, et al. Evaluation of impermeable covers for bedding in patients with allergic rhinitis. *N Engl J Med.* 2003 Jul 17; 349(3):237-46.
14. Frederick JM, Warner JO, Jessop WJ, Enander I, Warner JA. Effect of a bed covering system in children with asthma and house dust mite hypersensitivity. *Eur Respir J.* 1997 Feb; 10(2):361-366.
15. Woodcock A, Forster L, Matthews E, Martin J, Letley L, Vickers M, et al. Control of exposure to mite allergen and allergen-impermeable bed covers for adults with asthma. *N Engl J Med.* 2003 Jul 17; 349(3):225-236.
16. Cloosterman SG, Schermer TK, Bijl-Hoffland ID, Van Der Heide S, Branckeeff B, Van Den Elshout FJ, et al. Effects of house dust mite avoidance measures on Der p 1 concentrations and clinical condition of mild adult house dust mite-allergic asthmatic patients, using no inhaled steroids. *Clin Exp Allergy.* 1999 Oct; 29(10):1336-1346.
17. Dharmage S, Walters EH, Thien F, Bailey M, Raven J, Wharton C, et al. Encasement of bedding does not improve asthma in atopic adult asthmatics. *Int Arch Allergy Immunol.* 2006; 139(2):132-138.
18. Halken S, Hest A, Nikiassen U, Hansen LG, Nielsen F, Pedersen S, et al. Effect of mattress and pillow encasings on children with asthma and house dust mite allergy. *J Allergy Clin Immunol.* 2003 Jan; 111(1):169-176.
19. Luczynska C, Trethewell E, Smeton N, Burney P. A randomized controlled trial of mite allergen-impermeable bed covers in adult mite-sensitized asthmatics. *Clin Exp Allergy.* 2003 Dec; 33(12):1648-1653.
20. Reisman RE, Maurello PM, Davis GB, Georgitis JW, DeMasi JM. A double-blind study of the effectiveness of a high-efficiency particulate air (HEPA) filter in the treatment of patients with perennial allergic rhinitis and asthma. *J Allergy Clin Immunol.* 1990 Jan; 85(6):1050-7.
21. Wood RA, Johnson EF, Van Natta ML, Chen PH, Eggleston PA. A placebo-controlled trial of a HEPA air cleaner in the treatment of cat allergy. *Am J Respir Crit Care Med.* 1998 Jul; 158(1):115-20.
22. Antoniceilli L, Bilò MB, Pucci S, Schou C, Bonifazi F. Efficacy of an air-cleaning device equipped with a high efficiency particulate air filter in house dust mite respiratory allergy. *Allergy.* 1991 Nov; 46(8):594-600.

**BASELINE SYMPTOMS
UN-CORRELATED WITH EXPOSURE**




Figure 2: Mean Baseline symptoms (±SEM) for subjects by exposure classification based upon dust sample allergen concentrations (below threshold, sensitizing, elevated)⁴: (0-2,2-10,>10) ug/g for Mite group 1 and Can f 1; and (0-1,1-8,>8) ug/g for Fel d 1.

**TREATMENT EFFICACY
UN-CORRELATED WITH EXPOSURE**




Figure 3: Symptom reductions vs placebo (±SEM) for subjects by exposure classification based upon dust sample allergen concentrations (below threshold, sensitizing, elevated)⁴: (0-2,2-10,>10) ug/g for Mite group 1 and Can f 1; and (0-1,1-8,>8) ug/g for Fel d 1.

CLINICAL BENEFITS REPORTED FOR MULTIPLE OUTCOMES VS PLACEBO




Figure 4. Symptom Reductions and Quality of Life Improvements vs. Placebo: Percent change in total (TSS), nasal (TNSS), and ocular symptoms upon waking and quality of life (NRQLQ⁵ composite) from baseline. P-values indicate the significance of the Active – Placebo contrast.

CONCLUSIONS

The clinical benefits of PUREZONE were consistent across the study population and unrelated to dust sample allergen concentrations. This outcome suggests the following:

- Dust sample allergen concentrations may not have been representative of the actual allergen exposures faced by subjects in their bedrooms. This is possible because many patients benefited from night-time avoidance in bedrooms having allergen exposures classified as “below the threshold for sensitization.”
- Environments having allergen dust sample concentrations below the commonly accepted thresholds for sensitization⁴ may still provide enough exposure to cause symptoms. This is supported by baseline symptomatology, for which all subjects having bedroom exposures “below the threshold for sensitization” reported nocturnal symptoms, yet their only relevant sensitizations were to dog, cat, or dust mite allergen.

It is difficult to contrast our findings against those of other avoidance studies because:

- Trials of other personal air filtration devices did not involve allergen dust samples.
- Trials of other stand-alone avoidance measures (encasement, room HEPA) involving allergen dust samples have generally lacked clinical benefits. Multi-faceted trials typically involve targeting multiple allergen sources with multiple interventions that extend beyond the bedroom.